

NEW CortiCap Application Kits



Product Highlights

- KIT CONTAINS ALL SUPPLIES Ready to “grab and go”
- SINGLE-USE DISPOSABLE Eliminates concerns of cross-contamination
- AVAILABLE IN 3 SIZES Three sizes fit the wide range of head circumferences (Newborn, Child, and Adult)
- CONVENIENT AND COST-EFFECTIVE Package includes a QR Code link to Application Guide and Videos

The CortiCap™ is a disposable EEG electrode cap designed by Greentek, that provides the user an easy-to-use way to quickly apply scalp electrodes (positioned according to the International 10-20 System). The cap is comfortable for the patient to wear for extended recordings (48 hours or more).

The CortiCap Application Kit contains the supplies needed to prepare each electrode location and apply the conductive paste to produce high-quality EEG signals.





Designed by EEG Technologists for Today's Needs

The CortiCap was developed by EEG experts and produced in partnership with Greentek, a company known for quality electrodes and reusable electrode caps, who succeeded in designing a single-use cap without the compromises associated with earlier cap designs. The cap is made from a lightweight mesh material with silver/silver chloride ink printed electrodes on a flexible film attached to a silicone base, fixed on the cap. The cap covers the entire scalp and is held in place with a chin strap, securing the electrodes in position. The silicone base of each electrode cup is soft and shaped like a bowl with a wide skirt, providing comfort during EEG recording and is suitable for long-term monitoring. The patented silver/silver chloride ink printed electrodes guarantee high signal quality.



M-Quick Connector Cable
ORDERED SEPARATELY

Part Number	Description	Size	Head Circumference
KIT-101	CortiCap- Disposable EEG Electrode cap; Adult;	Adult	53 – 62 cm
KIT-102	CortiCap- Disposable EEG Electrode cap; Child	Child	45 – 53 cm
KIT-103	CortiCap – Disposable EEG Electrode Cap - Newborn	Newborn	28 – 36 cm
ACC-107	M-Quick Connector Cable for Adult/Child size		
ACC-108	M-Quick Connector Cable for Newborn size		

Ordering Information:

To order, go to www.corticare.com/corticap and link to the CortiCare eStore or download the CortiCap Order Form, complete and submit with payment (or Purchase Order with approved account). Contact us at estore@corticare.com for more information.

PACKAGE CONTENTS

Each Application Kit box contains:

CortiCap	Adult, Child, or Newborn
WavePrep	Small single-use cup
Cotton-tipped Applicators	6 each
Elefix conductive paste	180 g tube
20 mL "slip-tip" syringe	1 each

CortiCap Application Kits are shipped as a Case of 5 (five) boxes. An order of one case will contain 5 Application Kits.

The M-Quick Connector Cable is required to use the CortiCap and is ordered separately from the CortiCap Application Kits.

CortiCare, Inc.
5950 La Place Court, Suite 160
Carlsbad, CA 92008

estore@corticare.com

CortiCap is a trademark of CortiCare, Inc.
Greentek is a trademark of Wuhan Greentek Pty Ltd.